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**Directions of use**

# **P-Fill**

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Radeberg, den 6. 1. 2007

MEGADENTA Dentalprodukte GmbH

A handwritten signature in black ink, appearing to read 'Neumann', written in a cursive style.

Dr. M. Neumann

## Directions of use

### 1. Information

P-Fill Posterior-Filling Material

Shades: A1; A2; A3; A3,5; B2; D3

### 2. Composition

The inorganic filler particles comprise barium aluminium bore glass, spherical fillers, high-dispersal non-agglomerated silicon dioxide and agglomerated silicon dioxide.

The monomer matrix contains:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)-phenyl-propane,

7,9,9-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-di-yl-bis-methacrylate,

3,6-Dioxaoctamethylendimethacrylate,

Triethylenglycol-dimethacrylate.

### 3. Indications

For restorations of cavity classes I, II and V

### 4. Contraindication

Placement of P-Fill is contraindicated if the patient is known to be allergic to any of the ingredients in P-Fill.

### 5. Side Effects

Application of eugenol containing material, specifically base liners, is contra indicated, it will impair the polymerization.

### 6. Warning

The C-Cid etching contains 35% phosphoric acid. Avoid skin- and eyes contact with C-Cid and unpolymerized P-Fill material. P-Fill cures with light, it is sensitive to blue light (operating light, sun light). For that reason bright light during application of P-Fill should be avoided.

### 7. Processing and Placement

-Cavity Preparation

The cavity is prepared and dried in accordance with the general rules of composite filling treatment, the smear layer has not to be removed and unprepared dentine is cleaned with a cleaning paste. The prepared cavity is now rinsed with water and dried with oil free air.

-Base Liner

To avoid irritations resulting from direct contact to the dentin it should be shielded with a calcium hydroxide compound and a sub filling like glass ionomer cement.

-Cavity Finishing

The enamel margins are carefully finished with superfine diamond or carbide burs.

-Acid Etching

The cleaned enamel edges are etched with C-Cid. C-Cid is applied with the supplied brush.

After 20-40 s C-Cid is applied on the dentine allowed to take effect for another 10-15 s and then thoroughly rinsed with water. Excessive drying should be avoided in this case, as this would cause the exposed collagen fibre to collapse

-Application of the Primer

The Primer has now to be used after the recommendations of the producer.

Do not rinse off the dentin primer!

The cavity is dried with oil and water free air until the dentine surface appears dull and dry. It is essential that the primed dentine and the etched enamel are dry and contaminant free for the bond application.

C-Bond is applied onto the primed dentine in the dried cavity and the etched enamel, blown out thinly with oil free air and then cured for 10 s with a blue light source (e.g. Megalux).

Following the polymerization of C-Bond. Now P-Fill can be applied.

#### -Application

The quantity required for the cavity is extruded from the syringe and placed with an instruments into the cavity. The syringe should be capped and turned back immediately after extrusion. After modelling the occlusal surface the restoration has to be cured for 20-40 s with a blue light source (e.g. Megalux). The beam has to be directed to all surfaces of the restoration.

Restorations deeper than 2 mm should be filled in layers, each layer have to be cured separately (multi-phase technique). The influence of oxygen in the air causes a thin non-polymerized layer at the surface called "dispersion or smear layer". This layer should not to be removed as it will facilitate a solid chemical bond with the next layer of composite.

#### -Shaping and Finishing

The final forming and the removal of excess composite are carried out with rotating finishing instruments. For contouring and finishing fine and ultra fine finishing diamonds are used. For polishing flexible discs and polishing strips are suitable. All shaping and finishing have to be carried out under tooth preserving and with water-cooling. A sealing of the filling with C-Bond is recommendable.

### **8. Shelf life**

3 Years

### **9. Precautions and Storage**

Not to be stored above 25°C

Storage in refrigerator prolongs shelf life

### **10. Deliver units**

single syringe of 4,5g

part of assortments

### **11. Date information prepared**

January 2007

### **12. Manufacturer and Sales**

MEGADENTA Dentalprodukte GmbH

D-01454 Radeberg, Germany

---

1340 0107 516